

# 4 more Labor<sup>®</sup>

[WWW.4MORELABOR.COM](http://WWW.4MORELABOR.COM)



- **GLASS BEADS STERILIZERS GPS 150/300**

# GLASS BEADS STERILIZERS

## GPS 150/300

GPS 150 and GPS 300 are a fast, simple and accurate alternative to conventional sterilization methods.

The GPS 150 and GPS 300 glass bead sterilizers are equipped with a precision thermostat, a digital thermometer, a digital display and a temperature limiter. The sterilizers with a stainless steel chamber are compact enough to be used on Laminar Air Flow Work benches or other work tables.

In contrast to the conventional flame-removal method, they ensure a clean room atmosphere.

After switching on "ON", the devices are ready for operation within 25 minutes and ensure total destruction by sterilization.

All microorganisms in the sterilization chamber are killed within 10 seconds.

### FEATURES

	GPS 150	GPS 300
<b>Crucible size</b>	Φ 40 x 80	Φ 40 x 140
in mm		
<b>External dimensions</b>	130/145/155	130/145/215
Width / Depth / Height in mm		
<b>Temperature range</b>	100 bis 300	
in ° Celsius		
<b>Heating-up period</b>	≤25	
in Minuten		
<b>Weight of the glass beads</b>	150	300
in Gramm		
<b>Material of the lid</b>	High-grade steel	
<b>Order number</b>	10061869	10065540



# 4 more Labor<sup>®</sup>

WWW.4MORELABOR.COM | INFO@4MORELABOR.COM